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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

REF. 8104-8105-8106-8107 Knee orthosis for immobilization of 0°

REF. 8104F-8105F-8106F-8107F Knee orthosis for immobilization of 20°

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

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MAINTENANCE

-  Do not bleach
-  No chemical cleaning
-  Do not iron
-  Do not tumble-dry
-  Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- REF. 8104 – 8107:** Post-surgery or post-traumatic immobilization with knee in extension
- REF. 8104/F – 8107/F:** Post-surgery or post-traumatic immobilization with knee in extension of 20°

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Structure made of double-lined soft foam
- Terry cloth interior
- Medial and lateral rigid stays and back stays made of adjustable light alloy (only REF. 8104/F – 8107/F: with flexion of 20°)
- Adjustable medial and lateral panels
- Patella opening
- Washable

PUTTING ON THE APPLIANCE

ADAPTING FOR THE DOCTOR / TECHNICIAN

- Open the Velcro® tape closures and provisionally secure to each other. Loosen the side supporting plates from the rear plate.
- With the help of the supplied Velcro® tapes, place the rear supporting plate on the centre of the leg (fig. A). One end of the tape should be as close as possible to the edge of the plate.
- Secure the side supporting plates with the Velcro® closure to the rear plate and make sure that the rods are positioned precisely on the sides and the knee cap is located in the cut-out (fig. B). The opening on the front must lie precisely over the knee cap.
- If necessary, bend the rods to correspond to the anatomical profile of the leg.
- Tighten and secure the Velcro® closures using an alternating sequence.

PUTTING ON THE APPLIANCE FOR THE PATIENT

- Open the Velcro® tape closures and provisionally secure to each other.
- Close the Velcro® tapes on the rear supporting plate to provisionally secure to the knee immobilisation support. The opening on the front must lie precisely over the knee cap.
- Tighten and secure the Velcro® closures using an alternating sequence (high/low) starting with the closures closest to the centres of the knee (fig. C).

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